







# Characterisation

## Transmission Electron Microscopy

Karlsruhe Institute of Technology (KIT), Germany



<b>Contact:</b>	<b>Dr. Christian Kübel</b> Email christian.kuebel@kit.edu • Phone +49(721)608-28970 • Fax +49(721)608-28976						
<b>Material class:</b>	Silicon X	Polymer X	Metal X	Ceramic X	Glass X	Organic X	Other X
<b>Equipment:</b>	Titan 80-300 aberration corrected						
<b>Short technology description:</b>	<p>Transmission electron microscopy enables characterization of powders and thin films (which can be prepared in a target preparation from bulk materials) by direct imaging with up to atomic resolution. The image information can be locally correlated with spectroscopic techniques (EELS/EFTEM and EDX) to provide semi-quantitative elemental composition/maps with up to sub-nanometer resolution. All of these techniques can also be performed during in-situ heating or tensile testing of the material.</p> <p>For complex three-dimensional structures, electron tomography can be used to generate a 3D representation of the material with a spatial resolution of 1-2 nm, which can be used to quantitatively measure e.g. particle distributions, surface areas and faceting.</p>						
<b>Typical structures and designs:</b>			<p><b>BF-TEM image of Co/CoO core/shell particles with a diameter of ~30 nm</b></p>				
			<p><b>HRTEM image of Co/CoO particles showing the atomic arrangement in the single crystalline Co core and the polycrystalline CoO shell</b></p>				
			<p><b>EFTEM image of oxygen distribution in Co/CoO particles and local EELS spectrum of a CoO shell</b></p>				

		<p><b>Multivariate analysis of Co and O EFTEM maps revealing the Co und CoO<sub>x</sub> distribution</b></p>
		<p><b>Volume rendering of 3D tomographic reconstruction of Co/CoO particles</b></p>
		<p><b>Atomic resolution TEM image of a partial dislocation in nanocrystalline palladium</b></p>
<p><b>Special features:</b></p>	<ul style="list-style-type: none"> <li>- FEI Titan 80-300</li> <li>- (aberration corrected TEM)</li> <li>- Resolution:</li> <li>- 0.08 nm information limit TEM</li> <li>- 0.14 nm resolution in STEM</li> <li>- 0.7 eV energy resolution EELS</li> <li>- Techniques:</li> <li>- BF-TEM, aber. cor. HRTEM</li> <li>- HAADF-STEM, HRSTEM</li> <li>- EFTEM, EELS, EDX</li> <li>- Electron diffraction</li> <li>- Lorentz imaging</li> <li>- (S)TEM tomography</li> <li>- Low-dose techniques</li> <li>- In-situ heating holder (up to 900°C)</li> <li>- In-situ straining holder</li> </ul>	
<p><b>Limitations, constraints:</b></p>	<ul style="list-style-type: none"> <li>- Sample has to be a solid at RT and stable under vacuum conditions</li> <li>- Maximum sample thickness: 10-2000 nm (depending on resolution and technique)</li> <li>- Except in tomography, the TEM always provides an image/analysis of the projected structure of a sample</li> <li>- Depending on the chemical composition, the sample might be sensitive to the electron beam</li> <li>- H, He und Li can not be detected by analytical techniques</li> </ul>	
<p><b>Material examples:</b></p>	<ul style="list-style-type: none"> <li>- Thin films or nano powders can be directly imaged. All other materials need to be prepared for the TEM analysis</li> <li>- Target preparation by FIB lift-out (field of view typically 25x10 μm)</li> <li>- Classical preparation by cutting, grinding, argon ion milling or microtomy</li> </ul>	